

## OD600 Optimization Protocol

Some experiments require a specific cell density. Measuring cell density can be done by measuring the solution's absorption at 600nm wavelength with E.coli, as the culture solutions are roughly the same color.

Note that if the solution is too concentrated, math will be required to figure out how much distilled water to add. You probably won't survive.

### *Solutions*

Your favorite culture  
dH2O

### *Materials (sterile)*

Pipetman P20  
Pipetman P1000  
Spectrophotometry cuvettes

*You will also need access to a:*

Spectrophotometer

### *Procedure*

1. Swirl the **culture solution** flask gently to mix it and ensure uniform density.
2. Add **800µL culture** solution to a **cuvette** with a fresh, sterile **Pipetman P20 tip**.
3. Measure the OD600 of the cuvette with the **spectrophotometer**.
4. If the OD600 is too low, **incubate** the culture at **37° C** for **30 minutes** and take another reading. If the OD600 is too high, dilute with **XmL dH2O**, where:

$$X = \text{Current Volume} * \left( \frac{\text{Current OD 600}}{\text{Desired OD 600}} - 1 \right)$$

Use a **P1000** to do this unless the volume required is smaller than the minimum on the pipette.

Repeat this protocol until the OD600 is where you want it, making sure to mix the culture by swirling each time. Chill the cells on ice when the ideal optimization point is reached, as this will stymie their growth.