

# Pam's STET Plasmid Prep Protocol (Boiling Prep a la Birnboim and Doly)

*Modified 4/23/98; last modified 9/28/10*

## Reagents

### **STET Buffer**

8% sucrose  
50 mM TRIS-HCl, pH 8.0  
50 mM EDTA  
5% Triton X-100

Make up 100 ml STET Buffer at a time. *Autoclave*. NOTE: The Triton-X 100 often separates out after autoclaving; therefore, re-mix the phases well once solution has cooled down.

**10 mg/mL lysozyme** in STET buffer, FRESHLY-MADE – *Lately, I've been making a 20 mg/mL stock, but still adding 20 uL per sample (Step 4).*

### **Isopropanol**

### **70% ethanol**

### **TE**

## Protocol

1. Grow liquid cultures to saturation overnight at 37°C in 3-5 mL LB+Ampicillin or LB+Carbenicillin liquid (or patch out colony on 1/4 of LB+Amp or LB+Carb plate). Spin down ~1.5 mL of culture in a microfuge tube, 15 seconds at maximum speed. *If you want more DNA you can spin down a second aliquot of 1.5 mL in the same tube.*

2. Heat water bath to boiling (or make sure 100°C heating block is on).

3. Resuspend cells in 200 µL **STET buffer**, by gently pipetting up and down (or by “zotzing” and vortexing). *Lately, I've been adding 250 uL here.*

*NOTE: I put tubes on TOMY shaker to loosen pellets before adding STET buffer; in this case, pipeting up and down usually is unnecessary.*

4. Add 20 µL **lysozyme** solution. Mix briefly by quick vortex or inversion.

5. **Immediately** place at 100°C for 3 minutes.

*NOTE: If processing a lot of samples, it might be easier to resuspend loosened cell pellets directly in STET buffer containing lysozyme—with minimal vortexing—then immediately place in*

*boiling water bath.*

6. Remove promptly. Let tubes cool down for a few minutes on bench top. Spin tubes at top speed for 20 minutes.
7. Remove the gooey precipitate/pellet with a wooden toothpick and discard.
8. Add an equal volume of **100% isopropanol** to each tube (~250-300  $\mu$ L). Mix well. Put on ice or at 4°C for at least 30 minutes (RT is okay, too). If you are in a hurry this incubation can be shortened; however, the yield might be reduced. Also, this can be a stopping point.
9. Spin down precipitated nucleic acids (DNA+RNA) at top speed for 15 or 20 minutes. (The temperature is not important—pellets will be bigger at 4°C, but I think this is mostly due to salt.)
10. Decant or aspirate the supernatant and discard. Wash pellets with **70% ethanol**; vortex and respin for a few minutes.
11. Aspirate ethanol and dry pellet.

*NOTE: It is not necessary to dry the pellet, provided you aspirate all the ethanol.*

12. Resuspend pellet in 50-200  $\mu$ L **TE** (let sit at room temperature or on ice for a while).
13. Use 1-5  $\mu$ L per restriction digest. Don't forget to add RNase to digests or loading dye, as these preps will have RNA in them.
14. Store miniprep DNA at -20°C.