

Protocol: PCR Purification Protocol

QIAquick PCR Purification Kit (50)

All centrifugation steps are carried out at 13,000 rpm in a conventional table-top microcentrifuge at room temperature.

1. Add 5 volumes Buffer PB to 1 volume of PCR reaction and mix. If the color of the mixture is orange or violet, add 10 μ l 3 M sodium acetate, pH 5.0, and mix. The colour of the mixture will turn **yellow**.
2. Place a QIAquick column in a provided 2 ml collection tube.
3. To bind DNA, apply the sample to the QIAquick column and centrifuge for 1 min. Discard flow-through and place the QIAquick column back in the same tube.
4. To wash, add 0.75 ml buffer PE to the QIAquick column and centrifuge for 1 min. Discard flow-through and place the QIAquick column back in the same tube.
5. Centrifuge the QIAquick column once more in the provided 2 ml collection tube for 1 min to remove residual wash buffer.
6. Place each QIAquick column in a clean 1.5 ml microcentrifuge tube.
7. To elute DNA, add 50 μ l of distilled water (pH 7.0 to 8.5) or buffer EB to the centre of the QIAquick membrane and centrifuge the column for 1 min. For increased DNA concentration, add 30 μ l Buffer EB to the centre of the QIAquick membrane, let the column stand for 1 min, and then centrifuge.
8. If the purified DNA is to be analysed on a gel, add 1 volume of Loading Dye to 5 volumes of purified DNA. Mix the solution by pipetting up and down before loading the gel.

Buffers:

- Buffer PB - Binding buffer (has pH indicator added)
- Buffer PE - Wash buffer (has ethanol added)
- Buffer EB - Elution buffer

Note: Measure DNA concentration using NanoDrop Spectrophotometer before running any digestion