

Protocol : Generating Competent Cell

Before Starting

Pour minimal media plates.
 5x M9 salts etc overleaf.
 Prepare 100ml LB per strain.
 Prepare 50ml ice cold 0.1M CaCl₂ / 15% glycerol per strain
 Pre-chill eppendorf tubes

Day 1

1	Streak cells on minimal agar plate. Incubate 37°C overnight.
---	--

Day 2

2	Pick a colony into 5 ml LB + 100µl 1M MgSO ₄ Incubate 37°C in shaker overnight.
---	---

Day 3

3	Inoculate 100ml LB in pre-warmed conical with 1 ml of the 5 ml O/N culture from Day 2.
4	Incubate 2hrs in 37°C shaker until the cells at early log phase of growth curve ($A_{600} \sim 0.3$)
5	Transfer to chilled, sterile 50 ml Falcon tube and incubate on ice 10 min
6	Centrifuge at 3300g for 5 min at room temperature.
7	Resuspend in 10ml ice cold 0.1M CaCl ₂ / 15% glycerol and incubate on ice 30 min
8	Centrifuge at 3300g for 5 min at room temperature.
9	Resuspend in 1ml ice cold 0.1M CaCl ₂ / 15% glycerol. Transfer 100µl aliquots to pre-chilled, pre-labelled eppendorf tubes. Store at -70°C.

Solution Recipe

- 5x M9 salts in 500 ml dH₂O.

Na ₂ HPO ₄	32g	
KH ₂ PO ₄	7.5g	
NaCl		1.25g
NH ₄ Cl	2.5g	

- Minimal media plates for competence

In 50 ml Falcon;

39 ml	melted Bacteriological Agar solution ($\leq 50^{\circ}\text{C}$)
10 ml	5x M9 salt solution
1 ml	20% (w/v) D-glucose
50 μl	2mg ml ⁻¹ thiamine
5 μl	1M CaCl ₂
100 μl	1M MgSO ₄

- 0.1M CaCl₂ / 15% glycerol

In 50 ml Falcon;

5 ml	1M CaCl ₂
7.5 ml	100% glycerol