

Plasmid Mini-Preps using Qiagen solutions and No Columns

Marty Taylor – Boeke Lab. Last modified 12/2010. Based on Invitrogen Bacmid Isolation protocol and on Ed Davis's protocol from the Hughes lab. All spins at full speed in benchtop microcentrifuge at room temperature. It is important to line the lids up so all point the same way – that way you know where the pellets are!

1. Spin down 1.5/**2.0** ml of overnight bacterial culture or any cell suspension 30s in a microcentrifuge tube. Remove as much supernatant as possible (vacuum aspirator helpful)
 - a. Store glycerol stocks from leftover culture – I mix 666 μ L culture with 333 μ L 50% glycerol
2. Resuspend in 300/**400** μ l of Qiagen's P1 buffer containing RNase A by vortexing (or for fewer samples, pipetting works well). A multi-sample vortex head or Turbomix is very useful.
3. Add 300/**400** μ l of room temperature P2 (lysis) solution and mix by inversion. Incubate at room temperature for at least 5 minutes. The solution should change from opaque to almost clear. This incubation is critical for removal of RNA from the sample!
4. Add 300/**400** μ l of P3 (neutralization) buffer and mix by inversion (**do not vortex**). Buffer N3 can be used in place of P3 in a pinch, but the resulting DNA will contain more RNA and may work poorly with large plasmids. **Optional:** incubate on ice 5-10 min.
5. Centrifuge 5 min. If using N3 instead of P3, centrifuge 10-15min (pellet is less firm in N3).
6. **Optional:** Phenol/chloroform/Iso-Amyl alcohol extraction – cleans up the DNA, but expensive and we have found it doesn't make a difference for digests or sequencing, so I don't do it.
 - a. Transfer the supernatant to a new tube containing an equal volume of cold Ph/C/Iso-AA
 - b. Spin 5min at full speed at 4°C, take the top (aqueous phase) and proceed to IPA precip
7. Transfer 700-800/**1000** μ l supernatant into a new tube containing 800/**1000** μ l isopropanol. Avoid transferring the white pellet. Mix by inversion and incubate on ice at least 5 min
 - a. Pause point: can store at -20°C overnight or longer, but yield is not much changed
8. Centrifuge DNA precipitate 5 minutes. Discard supernatant; be careful – isopropanol pellets are only weakly adhered to the side of the tube and can be easily lost. Align the tubes before spinning!
 - a. Note: extending centrifuge time can increase yield
9. Rinse DNA pellet with 500 μ l of 70% EtOH. Centrifuge 5 minutes and aspirate all supernatant
 - a. Optional: repeat step 8.
 - b. Optional: centrifuge pellet again 10 seconds and remove supernatant
10. Air dry pellet at room temperature for 5-10 minutes.
11. Resuspend in 50 μ l of 10mM Tris, pH 8.0. TE is fine too, if the EDTA doesn't interfere with your next steps. Use 0.5 μ l in a 10 μ l restriction enzyme digest – typical yield is around 500-1500ng/uL in 50uL with high copy plasmids when starting with 1.5mL culture.

Buffer P1 (store at 4°C after RNase A added)

50 mM Tris-HCl pH 8.0 (6.06g/L)

10 mM EDTA (3.72g/L disodium salt dihydrate)

100 μ g/ml RNaseA (Qiagen 19051 or 19101 or we can make it)

Buffer P2:

200 mM NaOH (8 g/L)

1% SDS (10g/L or 50mL 20% solution)

Buffer P3: 3.0 M potassium acetate pH 5.5 (294.4g KOAc, ~110mL glacial acetic acid per L)