



UCL

Protocol: Transformation

1. Add 2 μ l of ice cold DNA to 50 μ l (1 vial) of competent cell on ice.
 2. Incubate on ice for 30 min.
 3. Heat shock cells by placing tubes in 42°C water bath for 45 seconds.
 4. Place cells on ice for 2 min.
 5. Add cells to 250 μ l LB medium (in a 15 ml or 50 ml tube)
 6. Shake for 1 hour at 37°C
 7. Place at least 50 μ l & 200 μ l of cells on separate selective nutrient agar plates. Spread 50 μ l plate first and then spread the 200 μ l plate.
 8. Let the cells grow overnight.
 9. Next day pick a colony into 2 ml selective LB. Can use 10 ml of LB at times to increase the cell sample.
-
9. Incubate at 37°C in shaker for 16 hours.

LB medium: lysogeny broth