

Transformation

Name:

Date:

Continue from
Experiment
Project Name:

Date

Name

Procedure

1. take cells from -80°C freezer and put them on ice! (every eppi contains about $400\ \mu\text{l}$ cells)
2. thaw cells on ice 20 minutes
3. pipette $50\ \mu\text{l}$ cells and $2\ \mu\text{l}$ DNA into eppi still on ice!
4. Incubate for 30 minutes on ice
5. Heat at 42°C for 60 sec
6. Incubate on ice for 5 minutes
7. Add $200\ \mu\text{l}$ LB Broth
8. Incubate for 2 hours at 37°C (cells with lysis cassette at $30^{\circ}\text{C}!!$)
9. Plate $50\ \mu\text{l}$ and $200\ \mu\text{l}$ on two different LB/Agar plates with appropriate antibiotic resistance

Documentation:

Why are you doing this experiment? Name of the sample? Where are they stored? Name the vector with inserts, antibiotika resistance etc.