

Ligation

Name:

Date:

Continue from

Date

Name

Experiment

Project Name:

Procedure

PCR tube:

total volume 20 μ l

1. add H₂O (17 μ l -X-Y-Z)
2. add 2 μ l Ligase Buffer 10x
3. add Insert 1, Insert 2 (when proceeding from 3A digestion use 2 μ l of each)
4. add Vector (20ng needed. When proceeding from 3A digestion use 2 μ l)
5. Add 1 μ l T4-DNA Ligase
6. Incubate 10-30 min at room temperature
7. heat for 20 minutes at 80°C
8. store at -20°C or directly proceed to transformation

Name of part	Ratio Insert:Vector = 3:1 or 1:1	Volume (μ l)
--------------	-------------------------------------	-------------------

X insert 1

Y insert 2

Z vector

H₂O

Documentation:

Why are you doing this experiment? Where are your parts stored? Name the parts for ligation etc.